

UCSSTM

Universal Cannulated Screw Set



- Cannulated instruments to aid placement of guide wires and screws
- Cannulated 4.0mm cortical screws are self-tapping for ease of insertion
- Threaded and non-threaded guide wires designed to facilitate screw placement
- Fully threaded and lag screw designs to maximize bone purchase
- MR and CT compatible



TITANIUM IMPLANTS (4.0mm Cannulated)


873-030 Cortical L=30mm	873-044 Cortical L=44mm	873-142 Cortical Lag* L=42mm
873-032 Cortical L=32mm	873-046 Cortical L=46mm	873-144 Cortical Lag* L=44mm
873-034 Cortical L=34mm	873-048 Cortical L=48mm	873-146 Cortical Lag* L=46mm
873-036 Cortical L=36mm	873-050 Cortical L=50mm	873-148 Cortical Lag* L=48mm
873-038 Cortical L=38mm	873-136 Cortical Lag* L=36mm	873-150 Cortical Lag* L=50mm
873-040 Cortical L=40mm	873-138 Cortical Lag* L=38mm	
873-042 Cortical L=42mm	873-140 Cortical Lag* L=40mm	*Partially threaded for fracture fixation


INSTRUMENTS

873-000	Implant/Instrument Case	873-007	Cannulated Drill Bit
873-001	UĈSS™ Guide	873-008	Cannulated Tap
873-002	Cannulated Trocar	873-009	Cannulated Screwdriver
873-003	Awl	873-010	T-Bar
000 004	0.00 0.	070 011	0411 Dt 11 0 t1 t /

873-005 Cannulated Depth Gauge

873-004 24" Disposable Guidewire (threaded) 873-011 24" Disposable Guidewire (non-threaded)





For further information and product availability, contact your MEDTRONIC SOFAMOR DANEK USA Sales Associate or call MEDTRONIC SOFAMOR DANEK USA Customer Service toll free: 800-933-2635.

Please refer to package insert for labeling limitations.

WARNING: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the



MEDTRONIC SOFAMOR DANEK USA FAX: 901-332-3920 1800 PYRAMID PLACE MEMPHIS, TN 38132 901-396-3133 www.sofamordanek.com

WATS: 800-876-3133 CUSTOMER SERVICE: 800-933-2635